


A Marker Selection analysis generates a heatmap that visualizes the differences in expression from one sample set to another.

Step 1: Select sample sets

To run a Marker Selection analysis, set up two sample sets in the Subset 1 and Subset 2 boxes on the **Comparison** tab. Subset 1 (1) is the reference group, and Subset 2 (2) is the comparison group. The example shows Subset 1 with LN Class V and Tub samples, and Subset 2 with LN Class III and Tub samples.

Click on the **Summary Statistics** (3) tab to review the number of samples that will be included in the analysis. Confirm that the number of samples available for each subset is sufficient for running a meaningful analysis.




Step 2: Set up analysis

Click on the **Advanced Workflow** tab (4) and select "Marker Selection (Differential Expression)" from the Analysis dropdown list (5).

Drag the high dimensional data node from the Navigate Concepts tree into the box (6). (Do NOT add a specific gene of interest, because the analysis will not run. If you would like to view a single probe, use the Heatmap analysis instead.)

Click the "Run" button (7). A popup box will show a progress bar as the data is compiled and analyzed by tranSMART. Some analyses, especially those with high dimensional data, may take several minutes to run.



Step 3: Review results

A heatmap and a table of top markers will appear underneath the variable selection box. The subsets are delineated by the orange and yellow bars at the top of the heatmap (8).

